



Your partner for
medical furniture



A Sense of Company

With more than 90 years of experience, SCHMITZ is one of the leading international manufacturers of medical products.

A family-owned business with a living tradition. Our typical employee has been with the company for some 18 years, proof of how strongly they identify with our values and philosophy.

In 2012, the company's management passed to Friedrich and Matthias Schmitz, the fourth generation of the family.


We are proud to be located in Germany and are deeply rooted in our home region. We develop and manufacture here.

90% of our suppliers are from Germany, of which 60% are from the local region.

Made in Germany



'The transition to the next generation is always something special.'



A Sense of Care

We understand the services which are our lifeblood, and aim to achieve maximum customer satisfaction. From consulting through planning to installation.


Cost-effectiveness through longevity is what you can expect from us.

With a well thought-out service concept and comprehensive support, SCHMITZ works alongside you as partner from start to finish.

Technical Customer Service

T +49 (0)2377 84 549

service@schmitz-medical.com



- Over 80 international trade partners
- Sales offices

SCHMITZ – represented worldwide







varimed®

**Couches, stools, multipurpose trolleys
and examination seats**

The indispensable helper in your daily routine.

Quality and longevity

High-quality materials combined with outstanding workmanship guarantee the longevity of our products.

The everyday helper

medi-net multipurpose trolley – our reliable helper that you can count on throughout your workday.







Examination couch 211.1080.0
upholstery gold orange 106, powder coating grey aluminium RAL 9007



ECG examination couch 211.1310.0
upholstery lavender purple 107, powder coating white aluminium RAL 9006



Examination and massage couch 211.9301.0
upholstery carmine red 92, powder coating grey aluminium RAL 9007



2 pcs. Side guards 200.4900.0




Examination couch 211.1180.0
upholstery ice blue 108, powder coating grey white RAL 9002

Examination couch 211.1180.0
upholstery agave green 109, powder coating grey white RAL 9002

varimed® Couches


¹ classic
² comfort

Model No. Description



211.1080.0 Examination couch¹²³

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 50 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°, 680 × 650 mm (H x W), safe working load: 225 kg.



211.1180.0 Examination couch¹²³

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°, 680 × 650 mm (H x W), safe working load: 225 kg.



200.0970.0 Paper roll holder

for paper roll up to 500 mm width



200.1030.0 Paper roll

500 mm wide

200.6500.0 Coloured frame

(w/o ill.) acc. to colour chart, for 211.1080.0 and 211.1180.0

Surcharge



200.4380.0 Head section with nose inlet

(w/o ill.) for 211.1180.0

Surcharge

200.4390.0 Filling for nose inlet


(w/o ill.)

Note:

For orthopaedic and massage treatment, couches with transverse and longitudinal bracings are recommended!

All examination couches with spacer block between frame and upholstery!

- ¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.
- ² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

varimed® **Couches**¹ classic
² comfort

Model No. Description

211.1310.0 ECG Examination couch¹²³ 

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with lateral bracing, 70 mm thick foam rubber upholstery with 2-section cover, levelling feet, head section adjustable up to +40°, 680 × 800 mm (H x W), safe working load: 225 kg.



200.0940.0 Paper roll holder
for 800 mm wide couches



200.1030.0 Paper roll
500 mm wide

200.5100.0 Longitudinal bracing
(w/o ill.) powder coated RAL 9002, RAL 9006 or RAL 9007,
for 211.1310.0

200.5120.0 Longitudinal bracing
(w/o ill.) powder coated according to colour chart

200.6500.0 Coloured frame
(w/o ill.) acc. to colour chart, for 211.1310.0

Surcharge

- ¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.
- ² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

Model No.	Description	
200.0000.0 (w/o ill.)	Different height up to max. 830 mm (from 700 mm height upwards, including transverse and longitudinal bracing) for couches 211.1180.0 – 211.1310.0	Surcharge
200.0010.0 (w/o ill.)	Mobility for couches 211.1180.0 - 211.1310.0, 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, longitudinal bracing. All castors are resistant to abrasions. (please indicate desired final height of couch!)	Surcharge
200.0020.0 (w/o ill.)	Different width⁴ of +/- 150 mm (max. width of 900 mm) for couches 211.1180.0 – 211.1310.0	Surcharge⁴
200.0030.0 (w/o ill.)	Different length⁴ of +/- 150 mm for couches 211.1180.0 – 211.1310.0	Surcharge⁴

⁴ If both special height and width are required,
the surcharge will only be applicable once.




Bariatric couch 211.1330.1
upholstery caribbean blue 82, frame chromium-plated

varimed® Couches

¹ classic
² comfort

Model No. Description



211.1330.1 Bariatric couch¹²³

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with transverse and longitudinal bracing, 2-section mattress, 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°, 680x800 mm (HxW). Safe working load: 300 kg. Mobile version and special widths/lengths not available!



200.0940.0 Paper roll holder for 800 mm wide couches



200.1030.0 Paper roll 500 mm wide


- ¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.
- ² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

varimed® 

Mobile couches, recovery couches

¹ classic
² comfort

Model No. Description



211.9301.0 Examination and massage couch¹²³

with transverse and longitudinal bracing, rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, 70 mm foam rubber upholstery with 2-section cover, head section adjustable up to +40°, 680 × 650 mm (H x W), safe working load: 225 kg.

All examination couches
with spacer block
between frame and
upholstery!

- ¹ 4 colours for leatherette upholstery in “classic” quality without surcharge.
- ² 15 colours for leatherette upholstery in “comfort” quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007.
Please indicate desired colour upon order placement.




Model No. Description

200.0970.0 **Paper roll holder**
for paper roll up to 500 mm width.




200.1030.0 **Paper roll**
500 mm wide.



200.4900.0 **1 Side guard and 1 side rail** **piece**
side guard laterally adjustable on side rail,
chromium-plated, with clamp, fitted to
rectangular steel tube frame
couches, may be lowered down to height of upholstery.

200.0020.0 **Different width of 800 mm**
(w/o ill.) for model 211.9301.0. **Surcharge**



200.5500.0 **Storage basket**
formed to nest onto transverse or longitudinal
bracing, basket 520 × 490 × 140 mm,
chromium-plated.

200.6500.0 **Coloured frame**
(w/o ill.) acc. to colour chart, for 211.9301.0. **Surcharge**

All castors are non-marking and resistant to abrasions.




varimed® C€^{MD}

Height adjustable examination and treatment couches

¹ classic
² comfort

Model No. Description



211.4100.0 Examination and treatment couch¹²³

hydraulic height adjustment 500 to 920 mm by means of foot lever. Servo-assisted head section adjustment from -20° to $+40^{\circ}$. 2-section upholstery 70mm thick, leatherette cover acc. to SCHMITZ colour chart. Powder coated chassis in white aluminium RAL 9006, grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover. With 4 levelling feet. Safe working load: 225 kg

211.4110.0 Examination and treatment couch¹²³

800 mm wide, otherwise like 211.4100.0.



211.4200.1 + 200.0250.0 + 200.0400.0

211.4200.1 Examination and treatment couch¹²³

electromotive height adjustment from 500 to 920 mm by means of foot control unit, with safety ON/OFF switch to turn control element on or off, 2.40 m connecting cable with mains plug, otherwise like 211.4100.0.

Upon request, hand control instead of foot control unit. For safety reasons, the height adjustment function is triggered by double-clicking on the hand or foot control.

211.4210.1 Examination and treatment couch¹²³

800 mm wide, otherwise like 211.4200.1.

- ¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.
- ² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

varimed® C€^{MD}

Height adjustable examination and treatment couches

Model No. Description



200.0250.0 Mobility

4 closed double castors Ø 100 mm,
2 of which el. conductive, central lock per axis,
1 castor with directional lock.
(Height adjustment range remains the same.)

Surcharge

**For mobility for retrofitting,
please contact our service team.**



200.0400.0 Side guards, pair

length: 990 mm, folding down to below top of
mattress, powder coated white aluminium RAL 9006.



200.0006.0 Holder for foot control

for Examination and Treatment couches 211.4200.1
and 211.4210.1, of steel, powder coated white
aluminium RAL 9006, upside-down mounting to
prevent accidental triggering

200.0930.0 Paper roll holder

(w/o ill.) for paper roll up to 590 mm wide.

200.1030.0 Paper roll

(w/o ill.) 500 mm wide

200.0410.0 Edge protection

(w/o ill.) of plastic, grey, at all four upholstery corners

200.0500.0 Head section with nose outlet

(w/o ill.)

Surcharge

200.0510.0 Filler for nose outlet

(w/o ill.)

Surcharge

Model No.	Description	
200.4402.0 (w/o ill.)	Side rail (piece) for Examination and Treatment couches 211.4100.0, 211.4110.0, 211.4200.1 and 211.4210.1, mounting either side, 25 × 10 × 750 mm, max. charge 30 kg, cannot be combined with side guards 200.0400.0. for 211.4100.0, 211.4110.0, 211.4200.1 and 211.4210.1.	
200.0251.0 (w/o ill.)	ECG shield with one potential equalization sockets each LHS and RHS as well as screening plate to reduce disturbance, suitable for 650 mm wide couches.	
200.0252.0 (w/o ill.)	ECG shield for 800 mm wide couches, otherwise like 200.0251.0.	
200.6530.0 (w/o ill.)	Surcharge for coloured chassis acc. to colour chart.	Surcharge



222.2760.0

222.2810.0

222.1550.0

222.1570.0

222.2760.0 upholstery milky coffee 102
222.2810.0 upholstery raspberry red 103
222.1550.0 and 222.1570.0 upholstery caribbean blue 82



221.1860.0
upholstery classic silver grey 55, powder coating grey white RAL 9002




242.2351.0

varimed® Swivel stools, Examination seats, Medical furniture

¹ classic
² comfort

Model No. Description




222.2810.0 Swivel stool¹²

5-arm plastic base, black, chromium-plated column, friction braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, chromium-plated column, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.

200.8000.0 Powder coated base according to colour chart
chromium-plated column.

Surcharge



222.2760.0 Swivel stool¹²



5-arm polished aluminium base, chromium-plated column, friction-braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.

¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.

² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.

varimed® Swivel stools, Examination seats, Medical furniture

¹ classic
² comfort



Model No. Description

222.1550.0 Examination stool¹²

specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release. 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection.

Safe working load: 135 kg.

222.1570.0 Examination stool¹²

specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release, foam rubber backrest with leatherette cover, 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection.

Safe working load: 135 kg.

200.0011.0 Brake loaded castors


braked when loaded, released when you come off the stool. For examination seats 222.1550.0 and 222.1570.0.

Surcharge

¹ 4 colours for leatherette upholstery in "classic" quality without surcharge.

² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.

¹ classic
² comfort




Model No. Description

221.1860.0 Examination/Blood sampling chair¹²³ €^{MD}
base of rectangular tube, powder coated in RAL 9002, RAL 9006 or RAL 9007, seat, backrest and armrests upholstered in leatherette classic black 57. Unassembled delivery. Safe working load: 225 kg.

200.6610.0 Special coloured finish of frame
acc. to colour chart.

Surcharge




242.2351.0 Drip-feed stand €^{MD}
5-arm base of 18/10 stainless steel frame, drip glass, double plastic castors diam. 50 mm, 2 of which lockable, antistatic, upper section with 4 18/10 stainless steel hooks, drip glass, height adjustment from 1,350 to 2,150 mm, one-hand action, footprint 635 mm. Safe working load: 2 kg (per hook).

- ¹ 3 further colours for leatherette upholstery in "classic" quality without surcharge.
- ² 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- ³ Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.



medi-net general purpose table 233.6060.0,
frame grey aluminium RAL 9007

varimed® medi-net general purpose table



Model No. Description

233.6060.0 medi-net general purpose table³
 rectangular steel tube frame, collapsible, powder coated grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, drawer with pre-assembled black galvanized ball bearing precision runners, table top, shelf and drawer front in pearl white 60 plastic veneers, handle for drawer front powder coated in grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, 4 double castors diam. 75 mm.
 The medi-net table 233.6060.0 can be fitted with 2 more drawers, in addition to the one standard drawer.

Unassembled delivery.

medi-net

Safe working load:

Totally: up to 50 kg

Drawer: up to 10 kg

Top/shelf: up to 20 kg

**201.3000.0 1 additional drawer
with handle**

200.6600.0 Coloured frame acc. to colour chart

Surcharge

200.6650.0 Coloured handle acc. to colour chart

Surcharge

³ Standard powder-coating colours for frame:
 9002, RAL 9006 or RAL 9007.
 Please indicate desired colour upon order placement.



medi-net general purpose table 233.6060.0
with accessories for injecting workplace:
201.1170.0, 201.1610.0, 201.1700.0, 201.2600.0



medi-net general purpose table 233.6060.0
with accessories for ergometry measuring
workstations 201.1300.0 (w/o equipment)



medi-net general purpose table 233.6060.0 +
201.0050.0 with accessories for home dialysis
201.1800.0, 201.1250.0 (w/o equipment)



Cannula and syringe dispenser 201.1610.0




Cannula and syringe dispenser 201.1650.0




varimed®

Accessories for medi-net general purpose table


Model No.	Description	
201.3000.0 (w/o ill.)	1 drawer with handle, running slides for subsequent fitting, (internal dimensions: 510 × 413 × 90 mm, W x D x H)	
201.0050.0	2 castors with locking device 75 mm Ø.	Surcharge
201.1300.0	Cable guide with fixing clamp 201.1910.0, 16 mm Ø.	
201.1170.0	1 round SST tray diam. 320 mm, with fixing clamp 201.1900.0, diam. 10 mm	
201.1250.0	Infusion pole CEM 2-section, with 4 hooks and fixing clamp 201.1920.0, diam. 18 mm, safe working load: 2 kg (per hook)	






Model No. Description


201.1610.0 Cannula and syringe dispenser
grey, first in/first out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0




201.1640.0 Cannula and syringe dispenser
grey, first in/first out system, without right twin set, requires fixing brackets 201.1700.0




201.1650.0 Cannula and syringe dispenser
grey, first in/last out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0



201.1660.0 Cannula and syringe dispenser
grey, first in/last out system, without twin set, requires fixing brackets 201.1700.0



201.1700.0 Fixing bracket, 1 pair
for use with 201.1610.0, 201.1640.0, 201.1650.0 and 201.1660.0



varimed®

Accessories for medi-net general purpose table

Model No. Description


201.1800.0 Console
SST 18/10 with fixing brackets

201.1900.0 Fixing clamp
(w/o ill.) 10 mm Ø


201.1910.0 Fixing clamp
(w/o ill.) 16 mm Ø

201.1920.0 Fixing clamp
(w/o ill.) 18 mm Ø

Model No. Description




201.2600.0 Ampule opener
"Cupfix"



201.2700.0 varimed Medical waste bin set
disposal set 1,5 l, 20 boxes made of polypropylene
with red lid, with holder

201.2710.0 varimed Medical waste bins, spare boxes
(w/o ill.) 1,5 l, carton of 100 boxes, without holder

201.2720.0 Holder
(w/o ill.) as 201.2700.0, however without boxes.
(Holders for different systems upon request)



201.2800.0 Refuse collector
with roll of bags, powder coated grey white

201.2850.0 Roll of spare bags
(w/o ill.) for 201.2800.0

varimed®

Model and material description

Couches

Tubular steel frame 40 × 20 mm, pair of feet 30 × 30 mm, powder-coated or chrome-plated. With levelling feet to compensate for floor unevenness. Head section with 2 adjustment mechanisms, 100% metal.

medi-net

Tubular steel frame 20 × 20 mm, plastic enamel coated.

Castors

All castors are non-marking and resistant to abrasions.

Drawers

Drawers run on precision ball-bearing slides.

Handles

of ergonomic design, distance of holes: 192 mm. The bowed bar handles are chemical resistant, easy to clean and to disinfect.

Powder coated finish

of high quality, multilayer gloss paint of extreme knock and scratch resistance. Easy to clean. Resistant to climatic conditions and disinfectants.

Standard finish of tubular steel frames

e.g. for couches: grey white RAL 9002, white aluminium RAL 9006, grey aluminium RAL 9007 or chromium-plated. Please refer to the different model descriptions.

Frames of couches, medi-net trolleys or stool bases

are also available in other colour acc. to colour chart at extra charge.

Further special powder coatings

acc. to RAL colour chart on request!

Upholstery

Foamed material on a support plate, with first-class, hard-wearing leatherette covers in two qualities.

Leatherette comfort (available for an additional charge): 15 colour tones

Wear-resisting material, with high resilience and flexibility providing excellent sitting and lying comfort. Grained surface.

Leatherette classic (available at no extra charge): 4 colour tones

Hard-wearing, flexible material. Smooth surface.

We recommend:**Accessories for upholstery cleaning**

Model No.	Description
2019878 (w/o ill.)	Special cleaner for leatherette 500 ml spray bottle (packaging unit 3 bottles)
2026973	Cleaning set for leatherette 1 × 500 ml spray bottle and a brush

Please refer to our separate operating manual for upholsteries for detailed information on cleaning and disinfection of upholsteries.

Leatherette comfort




gold yellow 86



agave green 109



gold orange 106



milky coffee 102



red carmine 92



sand beige 101



raspberry red 103



stone grey 96



ice blue 108



shale grey 97



caribbean blue 82




titan grey 91



atlantic blue 83



blossom white 84



lavender purple 107

Leatherette Powder coated finishes classic



white 53








silver grey 55



ocean 52



black 57

grey white
RAL 9002pure white
RAL 9010water lily white
RAL 0758510light saffron orange
RAL 0608530signal yellow
RAL 1003vermilion
RAL 2002lotus pink
RAL 3608510tricot lilac white
RAL 3108510tennis blue
RAL 2307020light blue
RAL 2606030sapphire blue
RAL 5003pallid blue
RAL 2308010lemon green
RAL 1008060filigree green
RAL 1207020ice-cold green
RAL 1209030white aluminium
RAL 9006grey aluminium
RAL 9007anthracite grey
RAL 7016

Plastic elements



silver-grey



antimicrobial

These colours have antimicrobial properties

Colours for upholsteries and plastic enamel coatings



Discover new, colourful worlds

With a total of 19 upholstery colours and 18 colours for powder-coated finish, you have many options to create colourful highlights in your surgery.

Whether muted or cheerful, classic or modern, elegant or extravagant – with our tool, you will certainly find the right colour for you in no time at all.

We would be delighted to advise you – please feel free to contact us.



Be inspired:
www.schmitz-medical.com



Contact

SCHMITZ medical GmbH
Zum Ostefeld 29
58739 Wickede (Ruhr)
Germany



All devices marked with the CE-MD symbol, in combination with our accessories, meet the requirements of Medical Device Regulation 2017/745/EU (MDR) as well as the Medical Device Amendment Act and are risk class I medical devices.



All devices marked with the simple CE symbol meet the requirements of low-voltage directive 2014/35/EU as well as the regulation on electrical equipment.

Please also refer to the instructions for use, for more details about the intended use!

T +49 (0)2377 84 0
F +49 (0)2377 84 162
zentrale@schmitz-medical.com

Exports department:
T +49 (0)2377 84 0
F +49 (0)2377 84 135
export@schmitz-medical.com

Technical service, hotline:
T +49 (0)2377 84 549
F +49 (0)2377 84 210
service@schmitz-medical.com

Information on
our establishments
worldwide can be found
on our website:

www.schmitz-medical.com



* This product is an article of merchandise, sold under OEM brand name.

Our goods and services are to be understood ex works. Deliveries are effected via our authorized partners or branch offices. The company reserves the right to alter the design, construction, dimensions and finish of its products without notice. All relevant dimensions and weights are to be understood as approximate information/data. Due to technical reasons colour deviations are possible. The articles shown in this catalogue are not authorized for sale in the United States.

Nominal charge 3.00 Euro.