



arco-matic® – The gynaecological examination and treatment chair

An award-winning design and state-of-the-art technology combined at excellent value for money.

The arco-matic® examination and treatment chair, with a width of upholstery of up to 680 mm at the back section and an individually adjustable head cushion with magnetic attachment, guarantees patients a pleasant, comfortable and safe sitting experience.

With an access height of only 550 mm, even older or physically challenged patients can easily take a seat on the chair, and even at its lowest level, the seat position remains horizontal.

Designed for a working load of 250 kg, the chair offers great safety and many treatment options even for adipose patients.

The ambient lighting and the new SCHMITZ colour concept combine to create a pleasant and anxiety-reducing atmosphere.

Configure your new arco-matic® now:
www.schmitz-medical.com



SCHMITZ – Your partner for gynaecology, video colposcopy, urology and proctology.

Made in Germany




Key Features



Horizontal position

The backrest can be moved to a horizontal position; this means that patients can be positioned flat, e.g. for an ultrasound examination, or in a shock position in case of an emergency.



Memory positions and pelvic elevation

An optimised adjustment speed and four freely selectable memory positions ensure quick adjustment to the relevant position. A lithotomy position is achieved by means of the pelvic elevation and various leg positioning systems.



Gynaecological cut-out and side rails

The deep gyn. cut-out of 155 mm offers improved access to the examination field. Additional side rails offer the possibility of attaching various accessories such as the LED examination light or various types of leg positioning systems.



Mobility

The mobility option turns the examination chair into a mobile unit – perfect for cleaning in the area under the arco-matic® and for flexible use in different treatment rooms.

A SENSE FOR COLOUR

SCHMITZ medical GmbH

Zum Ostfeld 29
58739 Wickede (Ruhr)

T +49 (0)2377 84 0

zentrale@schmitz-medical.com
www.schmitz-medical.com