



Cutting-edge technology combined with an award-winning design

The new medi-matic® 115.9 sets standards for use in urology!

The motorised back section, which can be raised up to 68°, not only provides a comfortable sitting position for the patient, but also optimises your urodynamic workflow.

Our gyro flushing system offers a direct water connection and drain and stays horizontal regardless of chair position. The extra drain gutter and the special seat section pad safely collect liquids so that they do not contaminate the chair column.

Closed surfaces without corners and edges ensure easy cleaning and disinfection. In addition, the mobility option turns the examination chair into a mobile unit – perfect for cleaning in the area under the medi-matic® and for flexible use in different treatment rooms.

A special highlight: the extremely low mounting height of only 470 mm and the maximum height of up to 1,175 mm.

SCHMITZ – Your partner for urology.



Configure your new medi-matic® now:
www.schmitz-medical.com

DISCOVER THE
medi-matic®

Application examples



Convenient examinations in urodynamics



Can be used as an outpatient surgical chair



Transurethral resection (TUR) procedure with Göpel-type leg supports



Benefit from the advantages of UniLeg leg supports




Rapid transition into a shock position




A removable rinsing basin, a special seat section pad and a seat section extension for optimised use in urodynamics


Highest position




Overview



Lowest position



Shock position



SCHMITZ medical GmbH

Zum Ostfeld 29
58739 Wickede (Ruhr)

Germany

T +49 (0)2377 84 0
zentrale@schmitz-medical.com
www.schmitz-medical.com